



Your local **gas generation** partner

Hydrogen Safety

A **PEAK**  **gas generation** brand

Why a generator is safe

- **Low stored gas volume** (<300cc) – Cylinders can contain >9000L
- **Low pressure operation** (up to 100psi) – Cylinders pressurized to 200+Barg
- **Real-time leak monitoring both internally and externally** – prevents excessive release of H₂ into lab or GC if a leak occurs
- **Days of supply following failure would be required to reach LEL**
– Cylinders contain sufficient gas volume to reach LEL in seconds

Consecutive failures required for critical build-up of H₂

1. **Hydrogen vent blockage** – causing H₂ build-up in generator
2. **Oxygen vent blockage** – causing O₂ build-up in generator
3. **Hydrogen exhaust blockage** – causing build-up of H₂ in PSA dryer
4. **Proportional valve failure** – preventing delivery of H₂ from generator
5. **Pressure sensor failure** – 3 sensors would need to fail simultaneously
6. **Valve failure** – 4 valves would need to fail simultaneously
7. **Non-return valve failure** – 3 non-return valves would need to fail
8. **Separation tank blockage** – Hydrogen forced into internal water tank

Ignition source would be required in addition to steps 1-8. System would also need to be hermetically sealed to reach the 4.1% LEL.

Single Cylinder failures required for Health risk

1. **Internal stress corrosion** cracking of the cylinder
2. **Cylinders can fail if over-pressurised** or weakened by the application of heat
3. **Cylinders are structurally weakened by the application of any heat source**
4. **Cylinders that may 'rocket'** if the regulator and valve assembly is damaged
5. **Cylinders protected by pressure relief valves**, fusible plugs or bursting discs that may not work correctly in a fire situation, or if damaged
6. **Accidental application of a compressed gas** or jet into eyes or into an open wound



**8 Consecutive
Failures
+ Combustion +
Sealed Environment**

VS



1 Failure

Consequences of a Cylinder Failure

The nature of the failure of a cylinder and its consequences depends on the combination of cylinder design, age and gas type. Flammable gases present greater risk but all cylinder failures can have significant safety consequences.

In a fire a cylinder can release its entire combustion energy. A cylinder of up to 1000L or 40kg bursting may result in a combination of the following:

- **Blast pressure** wave
- **Fireball** of up to 25 metres
- **A Projectile** traveling over 100 metres
- **Flying fragments** may travel up to 200 metres with high looping trajectories
- **Structural damage** to buildings in the vicinity

Generator safeguards

- **Full internal pressure decay test performed on start-up**
 - Any leak would prevent further production of H_2 and internal pressure would be vented.
Audible and visual alarms via HMI
- **Forced air ventilation**
 - Fans prevent build-up of H_2 and O_2 internally if leaks were to occur
- **Continuous monitoring of internal gas pressure**
 - Pressure sensor monitoring to prevent under/over pressure
- **Mechanical fail-safe in case of high pressure**
 - Automatic isolation and shutdown of the cell
- **When generator is in alarm status the following occurs:**
 - current to cell stopped
 - system depressurized
 - audible and visual alarms alert user of problem and fault information is displayed via the HMI



Gas generators provide a number of benefits



Convenience
Gas on-demand,
no cylinders to
change or
maintain supply
stocks



Consistency
Consistent gas
quality and supply,
no impurities or
running out of gas



Economy
Eliminate on-going
costs of cylinders,
manage lifetime
running costs



Safety
No pressurized
compressed gas
cylinders in your
lab




Green
No repeated gas
deliveries, energy
efficient

Hydrogen Generators A Safer Option for You, Your Colleagues, Your Workplace and the Environment

“there’s no need for us to have helium cylinders on site and no need to change air or helium cylinders at all, so that reduces instrument down time and cylinder manual handling risks.”


Ian Bennington, Senior Analyst, [Nerudia](#), UK

Safety - Hydrogen Cylinder vs Hydrogen Generator




A laboratory measuring 5m x 4m x 2.5m
has a volume of 50m³, or **50,000 litres**

The lower explosive level (LEL) of H₂ is only **4.1%**




We would only need **2050 litres of H₂** to reach the lower limit of LEL



A **50L gas cylinder** contains around **9000 litres** of hydrogen. Therefore releasing **25%** of the contents would reach **LEL**

A **Precision H₂ generator** produces up to **500 cc/min** and would take **almost 3 days** to reach the LEL (and this assumes no escape of hydrogen during this time)



Hydrogen Safety



“I recommend Peak Scientific generators to other laboratories due to safety issues, and eliminating the possible risks of using and handling cylinders.”

Askin Kayabasli, Laboratory Manager, Era Pharma, Turkey

Safeguard your hydrogen supply by using a generator

“a key benefit for us has been the elimination of the safety hazards associated with compressed gas cylinder handling.”

Kerri Heckrow, NPD Lab Manager, Evergreen Packaging, NC, USA

Certification



Evaluated and Certified to ensure compliance with both EMC and electrical safety of the system in the laboratory environment



Certified to ensure it can be safely operated, without risk of radiation frequencies causing interference/damage



Peak hydrogen generators meet the stringent demands of one of the World's most respected safety marks



Ensures that both the generator and its packaging are designed to minimise any damage in transit so that it will arrive at the lab in full working order



Peak H₂ generators meet IEC standards for laboratory use.

Dilution tests showed that the risk of explosion does not exist, even if all internal H₂ leaked and fans are not operational

More on Lab Gas Safety




H2 Safety

Blog




Lab Safety Focus

Blog



Changing carrier gas from He to H2

App note



Hydrogen Lab Safety

Poster

Thank you

www.peakscientific.com

Our Vision:

To exceed the expectation of our Customers, Colleagues and Suppliers

Our Values:

- **Respect** customer, supplier, colleague
- **Restless** constantly striving to improve
- **Freedom with responsibility**
- **Fun & passion** in everything we do

