

Discontinuation Notice

Sciex MS Tables

Scope

The purpose of this document is to serve notice of the discontinuation of a product, to define the level of service support and where applicable, advise of any new or alternative models that replace the discontinued product.

Product(s) Discontinued

The following product(s) are affected by this discontinuation notice:

Sciex Part Number	Peak Part Number	Model
1031379	10-1016-3	MS Table 1B
5019355	10-0002	MS Table 2B
1031380	10-1016-3N	MS Table 1N
5019354	10-0001	MS Table Size 2N (Int. Gas Generator)
5029543	10-1026	MS Table 2N (Hi-Flow Int. Gas Generator)

Effective Date

Peak and Sciex have established that, with immediate effect, the older-style table variants listed above are no longer required and should be removed from the supply chain in favour of current models; discontinuation of these products is therefore effective as of **13/10/2015**.

All discontinued part numbers have now been removed from the vendor price list. Peak Scientific will no longer accept orders for these table models.

Replacement Products

Replacements for each of the discontinued models are shown in the table below. Please contact your Peak Scientific account manager should the requirement arise.

Sciex Part Number	Peak part number	Model	Part Numbers replaced
5032309	10-9900	MS Table 1B	10-1016-3 10-0002
5032310	10-9901	MS Table 1N (Integrated Gas Generator)	10-1016-3N 10-0001
5032311	10-1028	MS Table 2N (Hi-Flow Int. Gas Generator)	10-1026

Service and Support

Part numbers which contain integral gas generators should continue to be serviced throughout the life of the product. The following consumable items should be purchased and service performed on an annual basis:

Part Number	Part Description
08-0317	Infinity 103* Series Annual Service Kit
Annual service kit contains:	
00-0045	Filter Element for BA filter (XA) blue
00-0046	Filter Element for BA filter (DAC) black
00-0039 x 2	RAC filter element x 2

Tables with integral gas generators (10-1016-3N, 10-0001 and 10-1026) will be fully supported and spares parts will be available until **October 2020**, at the earliest.

Part numbers 10-1016-3 and 10-0002 are non-maintenance products which have no service/support considerations; if spare panels are required for cosmetic reasons, please contact your Peak Scientific account manager.

It should be noted that due to reduced order quantities of associated service and spare parts, Peak Scientific reserves the right to increase part and service costs after the product discontinuation date.