

Dedicated **gas solutions** for **Agilent LC-MS**

With two decades of experience, Peak Scientific develop market-leading gas generators for labs. Our generators are engineered to deliver the consistent flow and purity you demand, at the push of a button. What separates us from our competitors is the service we provide. With certified Peak engineers stationed across the globe, our commitment is to keep your generator running day in, day out.

Why a Peak **gas generator**?



Convenient

No changing cylinders or running out during analysis



Consistent

Consistent quality and purity of gas, no risk of contamination



Economic

Eliminate on-going costs of cylinders, manage lifetime running costs



Efficient

Streamline workflow, faster sample throughput



Green

Reduce carbon footprint by eliminating recurring gas deliveries



Protected

Comprehensive on-site warranty & service contracts

Contact us today to discover more!

Web: www.peakscientific.com/agilent

Email: discover@peakscientific.com

Tel: +1-866-647-1649

Agilent **Solutions**

Peak Scientific has engineered various dedicated generators meeting the specifications of Agilent's LC-MS instruments.

Genius XE Series

- Up to 70L/min of contaminant-free nitrogen
- Can supply up to 2 LC-MS instruments (Genius XE 70 models)
- Next-generation high performance premium compressors, engineered exclusively for Genius XE
- Self contained, plug and play, 12 months service interval



Genius SQ 24

- Best value generator on the market
- Suitable for 1 x Agilent LC-MS
- Up to 24L/min of analytical grade nitrogen
- Compressor based solution, no need for an external air supply

Halo

- Suitable for 1 X Agilent MP-AES instrument
- Two independent outputs, one of up 25L/min of nitrogen and one of 36.5L/min of dry air
- Self-contained, plug and play, 12 months service interval
- Exceptional purity and purity stability



Solaris XE

- Up to 35 L/min of instrument-grade nitrogen
- Suitable for 1 x Agilent LC-MS/MS
- LC-MS/MS solution for labs with external air supply

Infinity XE 60

- Flow rates ranging from 10-520 L/min
- In-field expansion of nitrogen capacity to meet growing lab demands
- LC-MS/MS solution for labs with an external air supply

[**PEAK Protected**][™]

Peak Scientific gas generators come backed by a 12 month warranty (apart from Genius XE which has a 2 year warranty*). Beyond this period, you can ensure that your investment continues to be [Protected] by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

*Year 2 of warranty subject to generator being serviced at end of year 1 by a Peak-approved agent in accordance with fixed annual maintenance schedule. For full terms and conditions visit: www.peakscientific.com/warranty-statement