

Dedicated **gas solutions** for **Thermo Fisher Scientific**

With two decades of experience, Peak Scientific develop market-leading gas generators for labs. Our generators are engineered to deliver the consistent flow and purity you demand, at the push of a button. What separates us from our competitors is the service we provide. With certified Peak engineers stationed across the globe, our commitment is to keep your generator running day in, day out.

Why a Peak **gas generator**?



Convenient

No changing cylinders or running out during analysis



Consistent

Consistent quality and purity of gas, no risk of contamination



Economic

Eliminate on-going costs of cylinders, manage lifetime running costs



Efficient

Streamline workflow, faster sample throughput



Green

Reduce carbon footprint by eliminating recurring gas deliveries



Protected

Comprehensive on-site warranty & service contracts

Contact us today to discover more!

Web: www.peakscientific.com/thermoscientific

Tel: +1-866-647-1649 (North America)

Email: discover@peakscientific.com

+44 141-812-8100 (Europe)



Dedicated Solutions

Peak Scientific has engineered various dedicated generators meeting the specifications of Thermo Fisher Scientific's analytical instruments.

Genius XE Series

- Up to 70L/min of contaminant-free nitrogen
- Can supply up to 2 LC-MS/MS instruments (Genius XE 70 models)
- Next-generation high performance premium compressors, engineered exclusively for Genius XE
- Self contained, plug and play, 12 months service interval



Corona 1010A

- Dedicated solution for Thermo Fisher Scientific's Corona Veo and Vanquish instruments
- Tried and tested membrane technology
- Built-in air compressor for a complete gas solution

Solaris XE

- Up to 35 L/min of analytical grade nitrogen
- Suitable for 1 x Thermo Fisher Scientific LC-MS/MS
- LC-MS/MS solution for labs with external air supply



Other Thermo Solutions

Infinity XE 60

- Up to 520 L/min of analytical grade nitrogen
- Optional second outlet for Dry Air
- LC-MS/MS solution for labs with an external air supply

[PEAK Protected][™]

Peak Scientific gas generators come backed by a 12 month warranty. Beyond this period, you can ensure that your investment continues to be [Protected] by our comprehensive generator care cover. (Genius XE models have a 2 year warranty*)

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

*Year 2 of warranty subject to generator being serviced at end of year 1 by a Peak-approved agent in accordance with fixed annual maintenance schedule. For full terms and conditions visit: www.peakscientific.com/warranty-statement