

## Supplier Approval

This letter is to verify the Peak Scientific nitrogen generators listed below are fully verified or recommended based on requirement (not yet verified) to use with selected SCIEX MS systems.

Application	Built in compressor	No compressor
Echo® MS System	Genius XE SCI 2	Infinity 1032*
Triple Quad™ 7500 LC-MS/MS System	Genius XE SCI 2 MS Bench (G) SCI 2	Infinity 1032*
X500 Series QTOF TripleTOF® 6600 6500+ Series 6500 Series TripleTOF® 5600+ 5500 series TripleTOF® 4600 4500 Series API 4000 Series API 3200 Series Triple Quad 3500	Genius 1024 MS Bench (G) SCI 1	Infinity 1031
3200 MD	Genius XE 35	Solaris XE
4500 MD and Topaz® System	Genius 1024	Infinity 1031 Hi-Flow

\* Meet requirement, but not yet verified

12/18/2020

**X** Patty Sun

Patty Sun  
Senior Product Manager, MS Accessories  
Signed by: Sun, Patty