

Genius SQ 24

Nitrogen for Single Quad LC-MS



Your local *gas generation* partner

Description (제품 개요)


Genius SQ24는 Single Quad LC-MS시스템 전용 고순도 질소 가스 제네레이터 입니다. 작은사이즈로 실험실용 테이블 아래 설치가 가능합니다. 시중에 판매되는 모든 종류의 Single Quad LC-MS 사양에 맞춰 분당 최대 24L 까지 유량 조절이 가능합니다.

PEAK사의 20년 이상의 제품 설계 및 제조 기술을 바탕으로 제작된 Genius SQ24 는 연구자에게 최고의 분석 결과를 제공할 수 있게 고순도 질소를 안정적으로 분석기기에 공급합니다.



Applications

- Single Quad LC-MS



Key Features (제품 특징)

- Single Quad LC-MS용 고순도 질소가스 생산
- 유량조절가능 (최대 분당 24 L)
- 최대 116psi 까지 압력조절 가능 (최대압력시 분당 유량 22L)
- 보다 작아진 사이즈로 실험실용 테이블 아래 설치가 가능
- 컴프레서 내장형으로 외부 에어 공급 불필요
- 20년 이상의 제품 설계로 입증된 기술력을

제품 사양

	Genius SQ 24		
	100V	120V	230V
전압오프셋			
최대 유량	Up to 24 L/min		
최대 압력	Up to 116 psi (@ 22 l/min max)		
출력단자	1 x 1/4" BSPP		
적정 습도	80% Non-Condensing		
최대 고도	2000m		
프탈레이트	<0.01µm		
이상 물질	NONE		
적정온도	NONE		
정격 전압	5°C (41°F) to 35°C (95°F)		
발열량	50/60Hz, 11A	60Hz, 10A	50/60Hz, 6A
제품크기	<1,100VA	1,200VA	<1,380VA
적정온도 BTU	2934 @ 50Hz 3446 @ 60Hz	3685 @ 60Hz	5016 @ 50Hz 4231 @ 60Hz
제품크기 (H x W x D)	610 x 600 x 750 mm 24.0 x 23.6 x 29.5"		
실질무게	83 kg / 183 lbs		
선적 무게	104 kg / 229.3 lbs		

주문정보

파트 번호	3301661 (230v) 3301660 (120v) 3301650 (100v)
파트 번호	
유지보수계약	visit: www.peakscientific.com/ordering

† Noise level expressed as SPL (Sound Pressure Level) measured at 1m from source in a reverberant chamber

[PEAK Protected]™

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed by a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

피크사이언티픽 코리아
서울시 중구 퇴계로 100,
스테이트타워 남산 6층
(우) 04631

서비스 연락처: 최주남 과장
이메일: discover@peakscientific.com
전화번호: 010-2007-2722

www.peakscientific.com

Peak Scientific's Quality
Management System
conforms to:
ISO 9001:2015



FM 515973



Product Certifications

